

C INNER:

1. Test connections and set up:
  - Measure CINNER at CIMON on the front panel (5V/div.) and at AUX. OUTPUT on the Breakout Box (100mV/div) using an oscilloscope (100uS/div).
  - Triggers will be used that are generated by write commands to the Monitor Card.Insure that a working Monitor Card is installed in the crate.
  - CI/CO Range note: Low = 499/90.499k, High = 499/18.299k
  - IMPLANT Range note: Low = 1/43.2k, High = 1/7.15k. The output pulse is so small that IMPLANT MON will be measured at pins 7 & 11 of U90 which is prior to the divider circuit. U90 is located on the bottom side, middle section of the top portion of the board.
  - Triggers are generated by PUSHING the Trigger button on the Monitor Card, or by sending write commands to the Monitor Card, or run the LED PULSE routine

2. Power up settings:

CIDAC	10V
CI RANGE	0(low)
CI PULSE	0(disabled)
Delay Trig	0(disabled)

3. Test preparation settings:

readfile ibapacap.macro

Set CI DAC to 10V:

Set to high range:

Turn pulse enable on:

Turn Delay Trigger off:

*CinnerDAC(10)*

*ciRange(1)*

*ciPulse(1)*

*trigDelay(0)*

CIDAC	10V
CI RANGE	1(high)
CI PULSE	1(enabled)
Delay Trig	0(disabled)

4. C INNER test:

Input settings in Step 3 before starting test.

Send triggers to IBAPACAP:

press the manual button on the monitor card

CIDAC	10V
CI RANGE	1(high)
CI PULSE	1(enabled)
Delay Trig	0(disabled)


or

module X

*LEDpulse(20)*

Note: X is the slot that the monitor card is in, to make changes to the blip card change the module setting back to the slot number that the card is in.

Record output waveform:  
Expected output is a ramped squarewave with a flattop of 400uS, with an amplitude of: CIMON=10V, AUX.=300mV.



5 . Range test:  
Input settings in Step 3 before starting test.

Select low Range: *ciRange(0)*  
Send triggers to IBAPACAP: module X  
*LEDpulse(20)*

CIDAC	10V
CI RANGE	0(low)
CI PULSE	1(enabled)
Delay Trig	0(disabled)

*Note: X is the slot that the monitor card is in.*

Does the amplitude of either output change? \_\_\_\_\_  
Which one and what is new amplitude? \_\_\_\_\_

6. DAC test:  
Input settings in Step 3 before starting test.

Set DAC to various settings from 0 to 10V and record the measurements from the front panel and from the breakout box(BOB): *CinnerDAC(x)*

Send triggers to IBAPACAP: module X  
*LEDpulse(20)*  
*Note: X is the slot that the monitor card is in.*

CIDAC	1-5V	
CI RANGE	1(high)	
CI PULSE	1(enabled)	
Delay Trig	0(disabled)	
DAC	Ci Mon	BOB
0		
2		
5		
7		
10		

7 . PULSE DELAY test:  
Input settings in Step 3 before starting test.

Note: Adjustment of R43 will delay the pulse from 363 mSeconds to 36.7 Seconds.

CIDAC	10V
CI RANGE	1(high)
CI PULSE	1(enabled)
Delay Trig	1(enabled)

Turn Delay Trigger on: *trigDelay*

Adjust R43 so the voltage measured at junction of R44 and R43 is 5V.

Send a trigger by using the trigger button on the Monitor Card. Observe the output and time how long the before the pulse is observed. \_\_\_\_\_

Adjust R43 so the voltage measured at junction of R44 and R43 is 0.05V.

Send a trigger by using the trigger button on the Monitor Card. Observe the output and time how long the before the pulse is observed. \_\_\_\_\_

C OUTER:

1. Test connections and set up:
  - Measure COUTER at COMON on the front panel (5V/div.) and at AUX. OUTPUT on the Breakout Box (100mV/div) using an oscilloscope (100uS/div).
  - Triggers will be used that are generated by write commands to the Monitor Card. Insure that a working Monitor Card is installed in the crate.
  - CI/CO Range note: Low = 499/90.499k, High = 499/18.299k
  - Triggers are generated by PUSHING the Trigger button on the Monitor Card, by sending write commands to the Monitor Card, or run the LED PULSE routine(tricky way).

2. Power up settings:

CODAC	10V
CO RANGE	0(low)
CO PULSE	0(disabled)
Delay Trig	0(disabled)

3. Test preparation settings:

readfile ibapacap.macro

Set CO DAC to 10V:

Set to high range:

Turn pulse enable on:

Turn Delay Trigger off:

*CouterDAC(10)*

*coRange(1)*

*coPulse(1)*

*trigDelay(0)*

CODAC	10V
CO RANGE	1(high)
CO PULSE	1(enabled)
Delay Trig	0(disabled)

4. C OUTER test:

Input settings in Step 3 before starting test.

Send triggers to IBAPACAP: press the manual button on the monitor card

or

CODAC	10V
CO RANGE	1(high)
CO PULSE	1(enabled)
Delay Trig	0(disabled)


module X

*LEDpulse(20)*

Note: X is the slot that the monitor card is in.

Record output waveform:

Expected output is a ramped squarewave with a flattop of 400uS, with an amplitude of:COMON=10V, AUX.=300mV.



5. Range test:  
Input settings in Step 3 before starting test.

Select low Range: *coRange(0)*  
Send triggers to IBAPACAP: module X  
*LEDpulse(20)* Note: X is the slot that the monitor card is in.

CODAC	10V
CO RANGE	0(low)
CO PULSE	1(enabled)
Delay Trig	0(disabled)

Does the amplitude of either output change? \_\_\_\_\_  
Which one and what is new amplitude? \_\_\_\_\_

6. DAC test:  
Input settings in Step 3 before starting test.

Set DAC to various settings from 0 to 10V  
and record: *CouterDAC(x)*

CODAC	1-5V
CO RANGE	1(high)
CO PULSE	1(enabled)
Delay Trig	0(disabled)

Send triggers to IBAPACAP: module X  
*LEDpulse(20)*  
Note: X is the slot that the monitor card is in.

DAC	Co Mon	BOB
0		
2		
5		
7		
10		

IMPLANT SECTION:

1. Test connections and setup:
- IMPLANT RANGE note: The output pulse is so small that IMPLANT MON will be measured at pins 7 & 11 of U90 which is prior to the divider circuit. U90 is located on the bottom side, middle section of the top portion of the board. The voltage dividers are: Low = 1/43.2k, High = 1/7.15k.
  - Monitor IMPLANT output at IMPLANT MON. on front panel.

2. Power up settings:

IMPLANTDAC OUT	10V
HEAT DC	0(off)
IMPLANT PULSE ENA	0(disabled)
IMPLANT RANGE	0(low)

**BLIP TESTING PROCEDURE**      **CDMS**      **8 Oct. 1998**  
**IBAPACAP: BOARD #** \_\_\_\_\_ **STATUS** \_\_\_\_\_ **DATE** \_\_\_\_\_

### 3. Test preparation settings:

**readfie implant.macro**

Set IMPLANT DAC to 5V:	<b><i>ImplantDAC(5)</i></b>
Turn HEAT DC off:	<b><i>pHeat(0)</i></b>
Enable pulse:	<b><i>pPulse(1)</i></b>
Select high range:	<b><i>pRange(1)</i></b>

IMPLANTDAC OUT	5V
HEAT DC	0(off)
IMPLANT PULSE ENA	1(enabled)
IMPLANT RANGE	1(high)

#### 4 . DC HEAT test:

**Input settings in Step 3 before starting test.**

Turn HEAT DC on: *pHeat*

Record the measurements of Implant Mon on the front panel and of Pins 7 and 11 of U90.

Set IMPLANT DAC to different values from 0V to 10V and record the measurements.

Are the measurements equal in amplitude:

IMPLANTDAC OUT		5V
HEAT DC		1(on)
IMPLANT PULSE ENA		1(enabled)
IMPLANT RANGE		1(high)
DAC	Implant Mon.	Pin7
11		Pin
0		
1		
3		
5		
7		
10		

## 5. Send Pulses:

**Input settings in Step 3 before starting test.**

Disable HEAT DC:                   pHeat(0)  
Enable pulse:                       ***pPulse***  
Send triggers to IBAPACAP: module X

**LEDpulse(20)** Note: X is the slot that the monitor card is in.

IMPLANTDAC OUT	5V
HEAT DC	0(off)
IMPLANT PULSE ENA	1(enabled)
IMPLANT RANGE	1(high)

Are pulses observed at the measuring points? \_\_\_\_\_

6. **RANGE test:**

**Input settings in Step 3 before starting test.**

Turn HEAT DC on: *pHeat*

Measure the outputs of U90 pins 6, 7.

Are the outputs equal? \_\_\_\_\_ yes, condition OK)

Measure the outputs of U90 pins 10, 11.

Are the outputs equal?\_\_\_\_\_ (yes, condition OK)

```

IMPLANTDAC OUT 5V
HEAT DC          1(on)
IMPLANT PULSE ENA 1(enabled)
IMPLANT RANGE    1 & 0

```

Select low range: *pRange(0)*  
Measure pins 6 and 10 of U90.

Did the outputs on pins 6 & 10go away? \_\_\_\_\_ (yes, condition OK)

7. Continuity test:  
Use a DMM to measure continuity from the junction of R233 and R236 to pin 18 of P1.  
Is there continuity? \_\_\_\_\_
- Use a DMM to measure continuity from the junction of R235 and R237 to pin 34 of P1.  
Is there continuity? \_\_\_\_\_