MIPP TPC

David Eric Miller LLNL / Purdue University


Fermilab, 22 July 2005


TPC Cluster Fitting Why?

- Better cluster fitting improves dE/dX determinations, vertex finding, gain calibration and charge estimation.
- Currently we use Clusters Center of Mass.
- Appropriate fits are Gaussian (horizontal) and Gamma (Vertical).

Fitting Successes


- Gaussian distribution over the TPC Pads.
- Gamma distribution in time (height) over pads. (Amplifier shaping)







Fitting Failures


- A quick check for convoluted clusters is excessive width over pads.
- Same check for time buckets.







Gaussian Fit Parameters










How much does this correct?



Gamma Fit Parameters








How much does this correct?

(Peak resolution yet to be done)



What's Next?

- Gain needs to be folded into the error estimation.
- Studies' methods need to be reworked for general use.
 - Speed improvement.
 - General TPC Data Clusters are handled differently. Software Merge is Required.
 - Cluster pileups need to be de-convoluted.